MedicoDx™ Extraction Reagent Set for use with MedicoDx™ Streptococcal Grouping Kit
(for *in vitro* diagnostic use)

PRODUCT CODES **PL046C-7 (65 tests)**
 PL046C-18 (165 tests)

INTENDED USE

MedicoDx™ Extraction Reagent Set, used in combination with MedicoDx™ Streptococcal Grouping Kit (PL041C) provides a rapid method for the serological identification of groups A, B, C, D, F and G of the Lancefield groups of beta-hemolytic streptococci grown on agar plates.

SUMMARY AND EXPLANATION

Clinical, epidemiological and microbiological studies have conclusively shown that the diagnosis of streptococcal infections based on clinical symptoms always requires microbiological verification⁽⁴⁾. Beta-haemolytic streptococci are the most frequently isolated human pathogens among the representatives of the genus *Streptococcus*. Nearly all the beta-haemolytic streptococci possess specific carbohydrate antigens (streptococcal group antigens). Lancefield showed that these antigens can be extracted in soluble form and identified by precipitation reactions with homologous antisera. Different procedures for extraction of streptococcal antigens are currently in use^(1,2,6,7,10,11). The MedicoDx™ Streptococcal Grouping Kit is based on liberation of specific antigen from bacteria cell walls by modified nitrous acid extraction. The extracted antigen in conjunction with particle agglutination offers a rapid, sensitive and specific method for identification of streptococcal groups A, B, C, D, F and G from primary culture plates.

PRINCIPLE OF THE TEST

The MedicoDx™ streptococcal grouping method involves chemical extraction of group specific carbohydrate antigens using specially developed nitrous acid extraction reagents. The extraction reagents 1 and 2 provided in the set contain a chemical substance able to extract the streptococcal group specific antigens at room temperature. Extraction Reagent 3 contains a neutralizing solution. The neutralized extracts can be easily identified using blue particles sensitized with purified group specific rabbit immunoglobulins. These blue particles agglutinate strongly in the presence of homologous antigen and will not agglutinate when homologous antigen is absent.

REAGENTS / MATERIALS SUPPLIED

PL046C-7 / PL046C-18

Each set is sufficient for 65 / 165 extractions suitable for typing using MedicoDx™ Streptococcal reagents PL031C, PL032C, PL033C, PL034C, PL035C or PL036C. Materials are supplied ready for use.

Extraction Reagent 1 (PL047C-7 / PL047C-18): One dropper bottle containing 7 / 18 mL of extraction reagent 1 with 0.098% sodium azide as preservative.

Extraction Reagent 2 (PL048C-7 / PL048C-18): One dropper bottle containing 7 / 18 mL of extraction reagent 2.

Extraction Reagent 3 (PL049C-7 / PL049C-18): One dropper bottle containing 7 / 18 mL of extraction reagent 3 with 0.098% sodium azide as preservative.

The components are sold as a set of three vials.

PRECAUTIONS

1. Do not use reagents after expiry date shown on product label.
2. Some reagents contain sodium azide. Sodium azide can react explosively with copper or lead if allowed to accumulate. Although the amount of sodium azide in the reagents is minimal, large quantities of water should be used if reagents are flushed down the sink.
3. The extraction reagents contain a mildly caustic agent. In case of skin contact, immediately wash the area with soap and copious amounts of water. If the reagent comes into contact with an eye, flush with water for at least 15 minutes.
4. Safety precautions should be taken in handling, processing and discarding all clinical specimens as a pathogenic organism may be present.
5. The kit is intended for *in vitro* diagnostic use only.
6. The procedures, storage conditions, precautions and limitations specified in these directions must be adhered to in order to obtain valid test results.

7. Some reagents contain materials of animal origin and should be handled as a potential carrier and transmitter of disease.

STABILITY AND STORAGE

All set components should be stored at 2-8°C. **Do not freeze.** Reagents stored under these conditions will be stable until the expiry date shown on product label.

MATERIALS REQUIRED BUT NOT PROVIDED

- Blue Particle Suspensions for groups A (PL031C), B (PL032C), C (PL033C), D (PL034C), F (PL035C) and/or G (PL036C) streptococci.
- Polyvalent Positive Control "ABCDG" (PL040C)
- Disposable cards with 10 test circles (PL092C-48)
- Disposable mixing sticks (PL091C)
- Inoculating loops, Pasteur pipettes, borosilicate glass test tubes 12mm x 75mm, timer.

The following ATCC strains are recommended for use in quality control;

Organism	Lancefield Group	Reference
<i>Streptococcus pyogenes</i>	group A	ATCC 19615
<i>Streptococcus agalactiae</i>	group B	ATCC 12386
<i>Streptococcus dysgalactiae subsp. equisimilis</i>	group C	ATCC 12388
<i>Enterococcus faecalis</i>	group D	ATCC 19433
<i>Streptococcus sp. Type 2,</i>	group F	ATCC 12392
<i>Streptococcus dysgalactiae subsp. equisimilis</i>	group G	ATCC 12394

TEST PROTOCOL

Reagents must be used in conjunction with MedicoDx™ Streptococcal Grouping Kit (PL041C). Refer to package insert (PL041C).

**LIMITATION OF THE PROCEDURE
PERFORMANCE CHARACTERISTICS
REFERENCES**

Refer to kit package insert (PL041C)

<p>PL047C-7 PL047C-18</p>	<p align="center"> Warning</p> <p>Harmful if swallowed.</p> <p>Wash hands and exposed skin thoroughly after handling. Dispose of contents/container in accordance with local regulation.</p>
<p>PL048C-7 PL048C-18</p>	<p align="center"> Danger</p> <p>Causes severe skin burns and eye damage.</p> <p>Do not breathe vapour/mist. Wear eye protection/protective gloves/protective clothing. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents/container in accordance with local regulation.</p>

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