**INTENDED USE**

MedicoDx™ Streptococcal-Select Grouping Kit provides a rapid method for the serological identification of groups A, B, C, D, F and G of the Lancefield groups of beta-hemolytic streptococci grown on agar plates.

**SUMMARY AND EXPLANATION**

Clinical, epidemiological and microbiological studies have conclusively shown that the diagnosis of streptococcal infections based on clinical symptoms always requires microbiological verification (4). Beta-haemolytic streptococci are the most frequently isolated human pathogens among the representatives of the genus *Streptococcus*. Nearly all the beta-haemolytic streptococci possess specific carbohydrate antigens (streptococcal group antigens). Lancefield showed that these antigens can be extracted in soluble form and identified by precipitation reactions with homologous antisera. Different procedures for extraction of streptococcal antigens are currently in use (1, 2, 6, 7, 10, and 11). The MedicoDx™ Streptococcal-Select Grouping Kit is based on liberation of specific antigen from bacteria cell walls by modified nitrous acid extraction. The extracted antigen in conjunction with particle agglutination offers a rapid, sensitive and specific method for identification of streptococcal groups A, B, C, D, F and G from primary culture plates.

**PRINCIPLE OF THE TEST**

The MedicoDx™ streptococcal grouping method involves chemical extraction of group specific carbohydrate antigens using specially developed nitrous acid extraction reagents. The extraction reagents 1 and 2 provided in the kit contain chemical substances able to extract streptococcal group specific antigens at room temperature. Extraction Reagent 3 contains a neutralizing solution. The neutralized extracts can be easily identified using blue particles sensitized with purified group specific rabbit immunoglobulins. These blue particles agglutinate strongly in the presence of homologous antigen and will not agglutinate when homologous antigen is absent.

**REAGENTS**

Each kit is sufficient for 65 streptococcal grouping tests. Materials are supplied ready for use.

**MATERIALS PROVIDED**

Customers select any five vials of particle suspension to accompany the MedicoDx™ Streptococcal-Select Grouping Kit:

|                                      |            |
|--------------------------------------|------------|
| Blue Particle Suspension Group A     | PL031C     |
| Blue Particle Suspension Group B     | PL032C     |
| Blue Particle Suspension Group C     | PL033C     |
| Blue Particle Suspension Group D     | PL034C     |
| Blue Particle Suspension Group F     | PL035C     |
| Blue Particle Suspension Group G     | PL036C     |
|                                      |            |
| Extraction Reagent 1 (PL047C-7)      | } PL046C-7 |
| Extraction Reagent 2 (PL048C-7)      |            |
| Extraction Reagent 3 (PL049C-7)      |            |
| Polyvalent Positive Control 'ABCDGF' | PL040C     |
|                                      |            |
| Mixing Sticks                        | PL091C     |
| Reaction Cards                       | PL092C-48  |

**Particle Suspensions:** Each vial contains 3.5 mL of blue particles coated with purified rabbit antibodies to Group A, Group B, Group C, Group D, Group F or Group G streptococci. The blue particles are suspended in buffer pH 7.4 containing 0.098% sodium azide as preservative.

**Polyvalent Positive Control:** One dropper bottle containing 3 mL of **ready to use** polyvalent antigens extracted from inactivated streptococci of Lancefield groups A, B, C, D, F and G.

**Extraction Reagent 1:** One dropper bottle containing 7 mL of extraction reagent 1 with 0.098% sodium azide as preservative.

**Extraction Reagent 2:** One dropper bottle containing 7 mL of extraction reagent 2.

**Extraction Reagent 3:** One dropper bottle containing 7 mL of extraction reagent 3 with 0.098% sodium azide as preservative.

**ALL COMPONENTS OF THIS KIT ARE AVAILABLE SEPARATELY FOR PURCHASE.**

**MATERIALS REQUIRED BUT NOT PROVIDED**

- Inoculating loops
- Pasteur pipettes
- Test tubes 12 x 75 mm
- Timer

**PRECAUTIONS**

1. Do not use reagents after expiry date shown on product label.
2. Some reagents contain sodium azide. Sodium azide can react explosively with copper or lead if allowed to accumulate. Although the amount of sodium azide in the reagents is minimal, large quantities of water should be used if reagents are flushed down the sink.
3. The extraction reagents contain a mildly caustic agent. In case of skin contact, immediately wash the area with soap and copious amounts of water. If the reagent comes into contact with an eye, flush with water for at least 15 minutes.
4. Safety precautions should be taken in handling, processing and discarding all clinical specimens as a pathogenic organism may be present.
5. The kit is intended for *in vitro* diagnostic use only.
6. The procedures, storage conditions, precautions and limitations specified in these directions must be adhered to in order to obtain valid test results.
7. Some reagents contain materials of animal origin and should be handled as a potential carrier and transmitter of disease.

**STABILITY AND STORAGE**

All kit components should be stored at 2-8°C. **Do not freeze.** Reagents stored under these conditions will be stable until the expiry date shown on product label.

**SPECIMEN COLLECTION AND PREPARATION OF CULTURES**

For specific procedures regarding specimen collection and preparation of primary cultures refer to a standard microbiology textbook. In general, a fresh (18-24 hr.) gram positive beta-haemolytic (5 % sheep blood agar) isolate of streptococcal colonies is assumed. One to four large colonies should be sufficient for grouping; however if the colonies are minute, an increased number of colonies (loopful) should be used.

**TEST PROTOCOL**

- All components should be at room temperature (20-28°C) prior to use.
1. Re-suspend the test particles reagents by gently inverting the dropper bottle several times. Examine the dropper bottles to ensure that the particles are properly suspended before use. Do not use if the particles fail to re-suspend.
  2. Label one test tube for each isolate to be tested.
  3. Add 2 drops of Extraction Reagent 1 to each tube.
  4. Select 1-4 beta-haemolytic colonies using a disposable loop and suspend them in the Extraction Reagent 1. If colonies are minute, pick several well isolated colonies to be tested such that Extraction Reagent 1 solution becomes turbid. In all cases the streptococcal colonies should be picked from an area which contains the least amount of contamination.
  5. Add 2 drops of Extraction Reagent 2 to each tube.
  6. Mix the reaction by tapping the tube with a finger for 5-10 seconds.
  7. Add 2 drops of Extraction Reagent 3 to each tube. Mix the reaction as in **step 6.**
  8. Dispense one drop of each blue particle suspension onto separate circles on the test card labelled for each isolate being tested.
  9. Using a Pasteur pipette, for each test place one drop of extract beside each drop of the particle suspension.
  10. Mix the blue suspension and the extract with the sticks provided, using the complete area of the circle. A new stick should be used for each test circle.
  11. **Gently** rock the card allowing the mixture to flow slowly over the entire test ring area.
  12. At one minute, under normal lighting conditions, observe for agglutination.

**OPTIONAL TEST PROTOCOL:** For testing of specific groups and not the full panel, a ratio of 1:1:1 drop of each Extraction Reagent can be used. This ratio is sufficient for up to three tests.

**QUALITY CONTROL PROCEDURES**

Routine quality control procedures for each MedicoDx™ lot involve testing of the blue particle suspension and extraction reagents with each streptococci group A, B, C, D, F and G using ATCC strains or equivalent as listed in this section. The extract from these strains will agglutinate with the homologous particle suspension. The Polyvalent Positive Control is used to test the individual latex reagents and should be performed with each shipment and new kit lot number received.

The following ATCC strains are recommended for use in quality control;

| Organism   | Lancefield Group | Reference  |
|--|------------------|------------|
| <i>Streptococcus pyogenes</i>                        | group A          | ATCC 19615 |
| <i>Streptococcus agalactiae</i>                      | group B          | ATCC 12386 |
| <i>Streptococcus dysgalactiae subsp. equisimilis</i> | group C          | ATCC 12388 |
| <i>Enterococcus faecalis</i>                         | group D          | ATCC 19433 |
| <i>Streptococcus sp. Type 2,</i>                     | group F          | ATCC 12392 |
| <i>Streptococcus dysgalactiae subsp. equisimilis</i> | group G          | ATCC 12394 |

## QUALITY CONTROL PROCEDURES continued

In addition, the following procedure is recommended to check the performance of the reagents:

- As a test of absence of auto agglutination, the blue particle suspensions should not show agglutination with normal saline solution.

## INTERPRETATION OF RESULTS

**Positive results:** Rapid strong agglutination of the blue particles within one minute with one of the particle reagents indicates specific identification of the streptococcal isolate. A weak reaction with a single blue particle suspension should be repeated using a heavier inoculum. The repeated test is considered positive if agglutination occurs with only one of the blue particle suspensions. Figure 1 illustrates a suggested scheme for grouping streptococci.

**Negative results:** No agglutination of the blue particles. If traces of granulation are seen in the test circle the test should also be regarded as negative.

**Inconclusive results:** If weak clumping or non-specific reaction (stringiness) is present in the test circle after one minute, the test should be repeated using a fresh subculture. If the same result is seen after retesting, biochemical testing should be performed to identify the isolate.

**Non-specific results:** On a rare occasion you may see agglutination with more than one group. If this occurs, please check the purity of the culture used to perform the test. If it looks pure, repeat the test and confirm the identification of the isolate with biochemical testing.

## LIMITATION OF THE PROCEDURE

- False negative or false positive results can occur if insufficient amounts of culture or extraction reagents are used.
- The kit is intended for use in identification of beta-haemolytic streptococci only. If alpha or non-haemolytic streptococci are tested, the identification should be confirmed by biochemical tests (5, 9) (Refer to suggested scheme for grouping streptococci).
- False positive reactions have been known to occur with organisms from unrelated genera, e.g. *Escherichia coli*, *Klebsiella* or *Pseudomonas* (3, 8). These are likely to non-specifically agglutinate all particle reagents.
- Some strains of Group D streptococci have been found to cross react with Group G antisera; these strains may be confirmed as Group D by the bile-esculin test. Some strains of *Enterococcus faecium* and *Streptococcus bovis* might be difficult to be grouped.
- Listeria monocytogenes* may cross react with the Group B and/or G Streptococcal particle suspensions. The catalase test may be performed to distinguish between *Listeria*, which are catalase-positive, and streptococci, which are catalase-negative. Gram staining and motility testing may be performed as further aids to differentiation.
- Some strains of *Streptococcus milleri* (*Streptococcus anginosus*) typically non-haemolytic possess A, C, F or G antigens and can give positive reaction with Strep A, C, F or G latex reagents. Morphology on blood agar and biochemical testing should be used to identify these organisms.

## PERFORMANCE CHARACTERISTICS

### A. Cross - reactivity studies:

The MedicoDx™ Streptococcal Grouping Kit was tested for cross-reactivity using 33 ATCC reference strains. The kit successfully grouped all streptococci containing Lancefield groups A, B, C, D, F and G (N=16). No cross-reactivity was observed during the testing of other streptococcal strains (n=7) nor of other non-streptococcal organisms (n=10).

### B. Clinical performance studies:


1. The MedicoDx™ Streptococcal Grouping Kit was evaluated as part of a comparison of five commercially available streptococcal grouping kits. The study was performed at the Northern General Hospital in Sheffield, England. All of the kits were challenged with a panel of 302 beta-haemolytic streptococci composed of 64, 67, 44, 55, 56 and 4 strains of Lancefield groups A, B, C, D, G and F respectively. The results showed that 12 of the strains failed to group with any of the kits tested. Of the remaining 290 strains the MedicoDx™ Streptococcal Grouping Kit correctly identified 286 (98%). The authors concluded that the MedicoDx™ Streptococcal Grouping Kit proved to be both accurate and rapid, with sensitivity and specificity of 99% and 100% respectively. Furthermore, the average time to agglutination was substantially less than that achieved by three of the other four kits evaluated. Additional data is available upon request.

2. A second performance study was carried out at a Health Centre in Ontario, Canada. In this study, 111 primary cultures were included (110 tested, 1 inadequate). All the strains were originally grouped by Lancefield precipitation reactions. All group D were further biochemically confirmed using a BE (bile esculin) and PYR (pyrrolidonyl aminopeptidase) assay protocol. The primary cultures were tested in parallel using the MedicoDx™ Streptococcal Grouping Kit and an alternative grouping kit. In this study, the overall agreement between MedicoDx™ and Lancefield results occurred with 109 of 110 isolates tested (99%), while overall agreement between the alternative kit and Lancefield results occurred with 106 of 110 isolates tested (96.3%). The 110 primary isolates used in this study included 15 group A, 40 group B, 13 group C, 4 group D, 11 group F, 12 group G and 15 non-groupable strains.

## REFERENCES

- Ederer, G.M., Herrmann, M.M., Bruce, R. Matsen, J.M. and Chapman, S.S. (1972). Rapid Extraction Method with Pronase B for Grouping Beta-Haemolytic Streptococci. Appl. Microbiol., 23, 285.
- EL Kholy, A., Wannamaker, L.W. and Krause, R.M. (1974). Simplified Extraction Procedure for Serological Grouping of Beta-Hemolytic Streptococci. Appl. Microbiol., 28, 836.
- Elliott, S.D. and Tai, J.Y. (1978). The Type-Specific Polysaccharides of *Streptococcus suis*. J. Exp.Med.,148, 1699.
- Facklam, R.R. (1980). Streptococci and Aerococci, Ch. 8 in Manual of Clinical Microbiology, 3rd Ed., Edited by Lennette, E.H. Balows, A., Hausler, W.J., and Truant, J.P. American Society for Microbiology, Washington, D.C. page 88-110.
- Facklam R.R. (1977). Physiological Differentiation of Viridans Streptococci. J. Clin. Microbiol., 5, 184.
- Fuller, A.T. (1938). The Formamide Method for the Extraction of Polysaccharides from Haemolytic Streptococci. Brit. J. Exp. Path., 19, 130.
- Maxted, W.R. (1948). Preparation of Streptococcal Extracts for Lancefield Grouping. Lancet, ii, 255.
- Nowlan, S.S. and Deibel, R.H. (1967). Group Q Streptococci. I. Ecology, Serology, Physiology and Relationships to Established Enterococci. J. Bact., 94, 291.
- Petts, D.N. (1984). Early Detection of Streptococci in Swabs by Latex Agglutination Before Culture. J. Clin. Microbiol., 19, 432.
- Rantz, L.A. and Randall, E. (1955). Use of Autoclaved Extracts of Haemolytic Streptococci for Serological Grouping. Stanford Med. Bull., 13, 290.
- Watson, B.K., Moellering, R.C. and Kunz, L.J. (1975). Identification of Streptococci. Use of Lysozyme and Streptomyces albus filtrate in the Preparation of Extracts of Lancefield Grouping. J. Clin. Microbiol., 1, 274.

Figure 1 SUGGESTED SCHEME FOR GROUPING STREPTOCOCCI

\* Some strains of group D have been found to cross-react with group G antisera. [Harvey, C.L. and McIlmurray, M.B (1984) Eur.J.Clinical Microbiol, 10,641].

|  |  |
|--|--|
| <p><b>PL047C-7</b><br/><b>PL047C-18</b></p>  | <p> <b>Warning</b></p> <p>Harmful if swallowed.</p> <p>Wash hands and exposed skin thoroughly after handling.<br/>Dispose of contents/container in accordance with local regulation.</p>  |
| <p><b>PL048C-7</b><br/><b>PL048C-18</b></p>  | <p> <b>Danger</b></p> <p>Causes severe skin burns and eye damage.</p> <p>Do not breathe vapour/mist.<br/>Wear eye protection/protective gloves/protective clothing.<br/>IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.<br/>IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.<br/>IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.<br/>Dispose of contents/container in accordance with local regulation.</p> |
| <p><b>PL031C</b><br/><b>PL032C</b><br/><b>PL033C</b><br/><b>PL034C</b><br/><b>PL035C</b><br/><b>PL036C</b></p> | <p> <b>Danger</b></p> <p>May damage fertility or the unborn child.</p> <p>Obtain special instructions before use.<br/>Wear eye protection/protective gloves/protective clothing.<br/>If exposed or concerned: Get medical advice/attention.<br/>Dispose of contents/container in accordance with location regulation.</p>   |

|   |                  |   |                                      |
|---|------------------|---|--------------------------------------|
|  | = Use by         |  | = In vitro diagnostic medical device |
|  | = Lot number     |  | = Contains sufficient for <n> tests  |
|  | = Catalog number |  | = Temperature limitation             |
|  | = Manufacturer   |  | = Consult Instructions for use       |

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#### BIBLIOGRAPHIE

1. **Ederer, G.M., Herrmann, M.M., Bruce, R. Matsen, J.M. and Chapman, S.S.** (1972). Rapid Extraction Method with Pronase B for Grouping Beta-Haemolytic Streptococci. Appl. Microbiol., 23, 285.
2. **EL Kholi, A., Wannamaker, L.W. and Krause, R.M.** (1974). Simplified Extraction Procedure for Serological Grouping of Beta-Hemolytic Streptococci. Appl. Microbiol., 28, 836.
3. **Elliot, S.D. and Tai, J.Y.** (1978). The Type-Specific Polysaccharides of *Streptococcus suis*. J. Exp.Med.,148, 1699.
4. **Facklam, R.R.** (1980). Streptococci and Aerococci, Ch. 8 in Manual of Clinical Microbiology, 3. Auflage., Herausgegeben von Lennette, E.H. Balows, A., Hausler, W.J., und Truant, J.P. American Society for Microbiology, Washington, D.C. S. 88-110.
5. **Facklam R.R.** (1977). Physiological Differentiation of Viridans Streptococci. J. Clin. Microbiol., 5, 184.
6. **Fuller, A.T.** (1938). The Formamide Method for the Extraction of Polysaccharides from Haemolytic Streptococci. Brit. J. Exp. Path., 19, 130.
7. **Maxted, W.R.** (1948). Preparation of Streptokokken Extracts for Lancefield Grouping. Lancet, ii, 255.
8. **Nowlan, S.S. and Deibel, R.H.** (1967). Group Q Streptococci. I. Ecology, Serology, Physiology and Relationships to Established Enterococci. J. Bact., 94, 291.
9. **Petts, D.N.** (1984). Early Detection of Streptococci in Swabs by Latex Agglutination Before Culture. J. Clin. Microbiol., 19, 432.
10. **Rantz, L.A. and Randall, E.** (1955). Use of Autoclaved Extracts of Haemolytic Streptococci for Serological Grouping. Stanfor Med. Bull., 13, 290.
11. **Watson, B.K., Moellering, R.C. and Kunz, L.J.** (1975). Identification of Streptococci. Use of Lysozyme and *Streptomyces albus* filtrate in the Preparation of Extracts of Lancefield Grouping. J. Clin. Microbiol., 1, 274.

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